

Product datasheet for **HP203320**

USP9X Human qPCR Primer Pair (NM_001039590)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8239
Forward Sequence:	GTGTCAGTTCGTCTTGCTCAGC
Reverse Sequence:	GCTGTAACGACCCACATCCTGA
ACCN:	NM_001039590 , NM_001039590.1 , NM_001039590.2 , BC014792 , BC046205 , BC054861 , BC063645 , BC172429 , BE925861 , BX956755 , NM_001039590.3
UniProt ID:	Q93008
Synonyms:	DFFRX; FAF; FAM; MRX99; MRXS99F
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)