

Product datasheet for **HP203351**

C4orf29 (ABHD18) Human qPCR Primer Pair (NM_001039717)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80167
Forward Sequence:	CAGCTCTCTTGCACTGGCTAGA
Reverse Sequence:	CCAGATGCTGTGGACCAAGACA
ACCN:	<u>NM_001039717</u> , <u>NM_001039717.1</u> , <u>BC034253</u> , <u>BC034253.1</u> , <u>BC058847</u> , <u>BC058898</u> , <u>BC075793</u> , <u>BC111400</u> , <u>BC128142</u> , <u>BC128143</u> , <u>BX119453</u> , <u>BX500590</u> , <u>NM_001039717.2</u>
Synonyms:	C4orf29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)