

Product datasheet for **HP203484**

MOSPD3 Human qPCR Primer Pair (NM_001040097)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 64598 |
| Forward Sequence: | CACAGCACCTGCCAAATACACG |
| Reverse Sequence: | GTCCTGGACATCATAGTGGCTG |
| ACCN: | NM_001040097 , NM_001040097.1 , BC005042 , BC005042.1 , BC011653 , BC011653.2 , BC009482 , BC028475 , BG164079 , BG707362 , BQ923161 , NM_001040097.2 |
| UniProt ID: | O75425 |
| Synonyms: | CDS3; NET30 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)