

## Product datasheet for **HP203518**

### PAQR3 Human qPCR Primer Pair (NM\_001040202)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 152559   |
| Forward Sequence: | CACTGGGTTTGGCTCAATGGAG   |
| Reverse Sequence: | AGTACCGCTCTGGGACTTTGGA   |
| ACCN:             | <a href="#">NM_001040202</a> , <a href="#">NM_001040202.1</a> , <a href="#">BC029604</a> , <a href="#">BC031256</a> , <a href="#">BC047510</a> , <a href="#">BX648346</a> ,<br><a href="#">NM_001040202.2</a>  |
| UniProt ID:       | <a href="#">Q6TCH7</a>   |
| Synonyms:         | RKTG   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)