

Product datasheet for **HP203521**

DCAMKL2 (DCLK2) Human qPCR Primer Pair (NM_001040260)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	166614
Forward Sequence:	CCAGCAGCCAACCTTTCTACTCC
Reverse Sequence:	GTCAAGCTCAGGTCCACCGTTT
ACCN:	NM_001040260 , NM_001040260.1 , NM_001040260.2 , NM_001040260.3 , BC019831 , BC032726 , BC172430
UniProt ID:	Q8N568
Synonyms:	CL2; CLICK-II; CLICK2; CLIK2; DCAMKL2; DCDC3; DCDC3B; DCK2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)