

Product datasheet for **HP203691**

SDHAF1 Human qPCR Primer Pair (NM_001042631)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	644096
Forward Sequence:	TGCTGCGCATCGAGTACCTGTA
Reverse Sequence:	GTCGTGAGGCTGGGAACCCAC
ACCN:	<u>NM_001042631</u> , <u>NM_001042631.1</u> , <u>NM_001042631.2</u> , <u>BC007205</u> , <u>BC064605</u> , <u>BC069054</u> , <u>BC137516</u> , <u>BC137517</u> , <u>BM695501</u>
UniProt ID:	<u>A6NFY7</u>
Synonyms:	LYRM8; MC2DN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)