

Product datasheet for **HP203762**

CEMP1 Human qPCR Primer Pair (NM_001048212)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	752014
Forward Sequence:	CCATCCTATCTCTTTGGACCTGG
Reverse Sequence:	CCTTGCTTACAGGTGCTGTCTT
ACCN:	<u>NM_001048212</u> , <u>NM_001048212.1</u> , <u>NM_001048212.2</u> , <u>NM_001048212.3</u> , <u>BC146925</u> , <u>BC146926</u> , <u>BC152783</u> , <u>BC156602</u> , <u>BM681118</u>
UniProt ID:	<u>Q6PRD7</u>
Synonyms:	CP-23; CP23
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)