

Product datasheet for **HP203766**

SMIM15 Human qPCR Primer Pair (NM_001048249)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	643155
Forward Sequence:	GTGTGTGTAAGGGGCTGGAGTA
Reverse Sequence:	AGCAAATGGGATTCTTTAACAACAG
ACCN:	<u>NM_001048249</u> , <u>NM_001048249.1</u> , <u>NM_001048249.2</u> , <u>NM_001048249.3</u> , <u>BC014242</u> , <u>BC047489</u> , <u>BC066550</u> , <u>BC103711</u> , <u>BC146876</u> , <u>BC146881</u> , <u>BX538015</u> , <u>BX648459</u> , <u>NM_001048249.4</u>
UniProt ID:	<u>Q7Z3B0</u>
Synonyms:	C5orf43
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)