

Product datasheet for **HP203772**

Somatostatin Receptor 3 (SSTR3) Human qPCR Primer Pair (NM_001051)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6753
Forward Sequence:	ATGAGCACCTGCCACATGCAGT
Reverse Sequence:	CTTCACCACGATGAGCAGGTAG
ACCN:	<u>NM_001051</u> , <u>NM_001051.1</u> , <u>NM_001051.2</u> , <u>NM_001051.3</u> , <u>NM_001051.4</u> , <u>BC042068</u> , <u>BC096829</u> , <u>BF513812</u> , <u>BF514979</u> , <u>BF515539</u> , <u>BU584980</u> , <u>NM_001051.5</u>
UniProt ID:	<u>P32745</u>
Synonyms:	SS-3-R; SS3-R; SS3R; SSR-28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)