

Product datasheet for **HP203787**

TNFRSF1B Human qPCR Primer Pair (NM_001066)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7133
Forward Sequence:	CGTTCTCCAACACGACTTCATCC
Reverse Sequence:	ACGTGCAGACTGCATCCATGCT
ACCN:	<u>NM_001066</u> , <u>NM_001066.1</u> , <u>NM_001066.2</u> , <u>BC052977</u> , <u>BC011844</u> , <u>BC042167</u> , <u>BT019927</u> , <u>NM_001066.3</u>
UniProt ID:	<u>P20333</u>
Synonyms:	CD120b; p75; p75TNFR; TBPII; TNF-R-II; TNF-R75; TNFBR; TNFR1B; TNFR2; TNFR80
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)