

Product datasheet for **HP203790**

beta II Tubulin (TUBB2A) Human qPCR Primer Pair (NM_001069)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7280
Forward Sequence:	TTGGGAGGTCATCAGCGATGAG
Reverse Sequence:	AGGCTCCAGATCCACCAGGATG
ACCN:	<u>NM_001069</u> , <u>NM_001069.1</u> , <u>NM_001069.2</u> , <u>BC018780</u> , <u>BC018780.1</u> , <u>BC001194</u> , <u>BC127873</u> , <u>BC141924</u> , <u>BC141925</u> , <u>BM846690</u> , <u>BP381259</u> , <u>NM_001069.3</u>
UniProt ID:	<u>Q13885</u>
Synonyms:	CDCBM5; TUBB; TUBB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)