

Product datasheet for **HP203791**

gamma Tubulin (TUBG1) Human qPCR Primer Pair (NM_001070)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7283
Forward Sequence:	CACTCAAGAGGCTGACGCAGAA
Reverse Sequence:	GGTTGATCTGGGAGAAGGATGG
ACCN:	<u>NM_001070</u> , <u>NM_001070.1</u> , <u>NM_001070.2</u> , <u>NM_001070.3</u> , <u>NM_001070.4</u> , <u>BC000619</u> , <u>BC000619.2</u> , <u>BG167819</u> , <u>BT019931</u> , <u>NM_001070.5</u>
UniProt ID:	<u>P23258</u>
Synonyms:	CDCBM4; GCP-1; TUBG; TUBGCP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)