

Product datasheet for **HP203792**

Thymidylate Synthase (TYMS) Human qPCR Primer Pair (NM_001071)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7298
Forward Sequence:	GGTGTTTTGGAGGAGTTGCTGTG
Reverse Sequence:	GGAGAATCCCAGGCTGTCCAAA
ACCN:	<u>NM_001071</u> , <u>NM_001071.1</u> , <u>NM_001071.2</u> , <u>BC013919</u> , <u>BC013919.1</u> , <u>BC002567</u> , <u>BC018858</u> , <u>BC083512</u> , <u>BT006811</u> , <u>NM_001071.4</u>
UniProt ID:	<u>P04818</u>
Synonyms:	HST422; TMS; TS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)