

## Product datasheet for **HP203802**

### TMUB2 Human qPCR Primer Pair (NM\_001076674)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79089
Forward Sequence:	CCTCATCACTGTGCGGCTCAA
Reverse Sequence:	CAGTTTCATCTGGCTTTCTGTCC
ACCN:	<a href="#">NM_001076674</a> , <a href="#">NM_001076674.1</a> , <a href="#">BC063489</a> , <a href="#">BC063489.1</a> , <a href="#">BC001792</a> , <a href="#">BC004306</a> , <a href="#">BI906577</a> , <a href="#">BM999307</a> , <a href="#">NM_001076674.3</a>
UniProt ID:	<a href="#">Q71RG4</a>
Synonyms:	FP2653
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)