

## Product datasheet for **HP203851**

### TSEN34 Human qPCR Primer Pair (NM\_001077446)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 79042  |
| Forward Sequence: | GGAGGTGACTTCCTGGTCTATC   |
| Reverse Sequence: | CTTTCTGACGCTGGTTCCAAGG   |
| ACCN:             | <a href="#">NM_001077446</a> , <a href="#">NM_001077446.1</a> , <a href="#">NM_001077446.2</a> , <a href="#">BC004530</a> , <a href="#">BC004530.2</a> , <a href="#">BC020805</a> , <a href="#">BC020805.1</a> , <a href="#">BC000944</a> , <a href="#">BE779123</a> , <a href="#">BM977824</a>  |
| UniProt ID:       | <a href="#">Q9BSV6</a>   |
| Synonyms:         | LENG5; PCH2C; SEN34; SEN34L  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)