

## Product datasheet for **HP203862**

### TCF19 Human qPCR Primer Pair (NM\_001077511)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 6941   |
| Forward Sequence: | TCAGCCTGGAAGACCACAGCAG   |
| Reverse Sequence: | CCAAAGGTCAGGAGGTCTCCAT   |
| ACCN:             | <a href="#">NM_001077511</a> , <a href="#">NM_001077511.1</a> , <a href="#">BC033086</a> , <a href="#">BC033086.1</a> , <a href="#">BC002493</a> , <a href="#">BC044632</a> , <a href="#">BP365938</a> , <a href="#">BT007247</a> , <a href="#">BX647174</a>   |
| UniProt ID:       | <a href="#">Q9Y242</a>   |
| Synonyms:         | SC1; TCF-19  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)