

Product datasheet for **HP203904**

ASAH2 (ASAH2B) Human qPCR Primer Pair (NM_001079516)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	653308
Forward Sequence:	GAGGCTACTTCAACATCGTGCC
Reverse Sequence:	TAGATTCCAGGCTGGCAGTGT
ACCN:	<u>NM_001079516</u> , <u>NM_001079516.1</u> , <u>NM_001079516.2</u> , <u>NM_001079516.3</u> , <u>BG773189</u> , <u>BI459539</u> , <u>NM_001079516.4</u>
UniProt ID:	<u>P0C7U1</u>
Synonyms:	ASAH2C; ASAH2L; bA98I6.3; bA449O16.3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)