

Product datasheet for **HP203925**

Pea3 (ETV4) Human qPCR Primer Pair (NM_001079675)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2118
Forward Sequence:	AGGAACAGACGGACTTCGCCTA
Reverse Sequence:	CTGGGAATGGTCGAGAGTTT
ACCN:	<u>NM_001079675</u> , <u>NM_001079675.1</u> , <u>NM_001079675.2</u> , <u>BC016623</u> , <u>BC016623.1</u> , <u>BC007242</u> , <u>BF304101</u> , <u>BM013535</u> , <u>BP268896</u> , <u>BU166857</u> , <u>NM_001079675.5</u>
UniProt ID:	<u>P43268</u>
Synonyms:	E1A-F; E1AF; PEA3; PEAS3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)