

Product datasheet for **HP204012**

KIAA1543 (CAMSAP3) Human qPCR Primer Pair (NM_001080429)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57662
Forward Sequence:	CACCATCCACTGCTATTGTCCC
Reverse Sequence:	AATCCTGCACGAGCTGGAGTT
ACCN:	<u>NM_001080429</u> , <u>NM_001080429.1</u> , <u>NM_001080429.2</u> , <u>BC020431</u> , <u>BC035808</u> , <u>BC036935</u> , <u>BC064490</u> , <u>BC128201</u> , <u>BC128202</u> , <u>BU733578</u> , <u>NM_001080429.3</u>
UniProt ID:	<u>Q9P1Y5</u>
Synonyms:	KIAA1543; NEZHA; PPP1R80
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)