

Product datasheet for **HP204024**

TTC36 Human qPCR Primer Pair (NM_001080441)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	143941
Forward Sequence:	GCTTCAGCCTACAACAACCGTG
Reverse Sequence:	CAGGAGTCCGCGCTGCACAAA
ACCN:	<u>NM_001080441</u> , <u>NM_001080441.1</u> , <u>NM_001080441.2</u> , <u>BC146916</u> , <u>BC146921</u> , <u>BC171906</u> , <u>BC171913</u> , <u>NM_001080441.4</u>
UniProt ID:	<u>A6NLP5</u>
Synonyms:	HBP21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)