

Product datasheet for **HP204036**

INTS1 Human qPCR Primer Pair (NM_001080453)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26173
Forward Sequence:	CTCTGTTGTCCATCTTCTCACGC
Reverse Sequence:	TGCTCCAGCGTAGGAAGACCTG
ACCN:	NM_001080453 , NM_001080453.1 , NM_001080453.2 , BC004286 , BC013367 , BC018777 , BC069262 , BF346658 , BK005720 , BM547069 , BM553299 , BX324321 , NM_001080453.3
UniProt ID:	Q8N201
Synonyms:	INT1; NDCAGF; NET28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)