

Product datasheet for **HP204047**

ASPG Human qPCR Primer Pair (NM_001080464)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	374569
Forward Sequence:	CTGTGACCACAGACTATGCAGC
Reverse Sequence:	AGCAGCTCCTTCCTGACATCCA
ACCN:	<u>NM_001080464</u> , <u>NM_001080464.1</u> , <u>NM_001080464.2</u> , <u>BC035836</u> , <u>BC136637</u> , <u>BM544197</u> , <u>BX247999</u> , <u>BX280385</u> , <u>BX407951</u>
UniProt ID:	<u>Q86U10</u>
Synonyms:	C14orf76; GPA/WT; LYSOLP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)