

## Product datasheet for **HP204056**

### MFSD2B Human qPCR Primer Pair (NM\_001080473)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388931
Forward Sequence:	ATCTGGCGTGAGCATTGCTGTG
Reverse Sequence:	ATGCGCCAGACAGCTTGGTGAA
ACCN:	<a href="#">NM_001080473</a> , <a href="#">NM_001080473.1</a> , <a href="#">BC033385</a> , <a href="#">BX357394</a> , <a href="#">BX369190</a> , <a href="#">BX647146</a>
UniProt ID:	<a href="#">A6NFX1</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)