

Product datasheet for **HP204063**

MBOAT1 Human qPCR Primer Pair (NM_001080480)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	154141
Forward Sequence:	GGTTTCCACAGCTTGCCAGAAC
Reverse Sequence:	ACCAGTCATCCACAAGGCAGGT
ACCN:	<u>NM_001080480</u> , <u>NM_001080480.2</u> , <u>BC045695</u> , <u>BC150650</u> , <u>BC150652</u> , <u>BC156593</u> , <u>BM978819</u> , <u>NM_001080480.3</u>
UniProt ID:	<u>Q6ZNC8</u>
Synonyms:	dj434O11.1; LPEAT1; LPLAT; LPLAT 1; LPSAT; OACT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)