

Product datasheet for **HP204105**

CC2D2A Human qPCR Primer Pair (NM_001080522)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57545
Forward Sequence:	GGAAGCGTAACACCCAATGACC
Reverse Sequence:	AGTGAAC TCGGAAGTCTGCTCC
ACCN:	<u>NM_001080522</u> , <u>NM_001080522.1</u> , <u>NM_001080522.2</u> , <u>BC053865</u> , <u>BC061921</u> , <u>BC070395</u> , <u>BC103710</u> , <u>BX647334</u>
UniProt ID:	<u>Q9P2K1</u>
Synonyms:	COACH2; JBTS9; MKS6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)