

Product datasheet for **HP204136**

CCDC32 Human qPCR Primer Pair (NM_001080791)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90416
Forward Sequence:	CGGTTCTCCAGGAGAAGTTAG
Reverse Sequence:	CCTGTGACTCTGGAGGAATCAG
ACCN:	<u>NM_001080791</u> , <u>NM_001080791.1</u> , <u>NM_001080791.2</u> , <u>BC001673</u> , <u>BC012189</u> , <u>BG709516</u> , <u>BP202606</u> , <u>BP289634</u> , <u>BQ878848</u> , <u>BU657334</u> , <u>BX404085</u> , <u>NM_001080791.3</u>
UniProt ID:	<u>Q9BV29</u>
Synonyms:	C15orf57
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)