

## Product datasheet for **HP204259**

### MZT2A Human qPCR Primer Pair (NM\_001085365)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	653784
Forward Sequence:	CGCAAGAAGGTGCTGAGCACC
Reverse Sequence:	GCTTCAGCAGGTCCACCAGGA
ACCN:	<u><a href="#">BC018206</a></u> , <u><a href="#">NM_001085365</a></u> , <u><a href="#">NM_001085365.1</a></u> , <u><a href="#">BC018206.1</a></u> , <u><a href="#">BC063024</a></u> , <u><a href="#">BC104650</a></u> , <u><a href="#">BI599271</a></u> , <u><a href="#">BM554057</a></u> , <u><a href="#">BU682101</a></u> , <u><a href="#">NM_001085365.2</a></u>
UniProt ID:	<u><a href="#">Q6P582</a></u>
Synonyms:	FAM128A; MOZART2A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)