

Product datasheet for **HP204282**

C17orf67 Human qPCR Primer Pair (NM_001085430)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	339210
Forward Sequence:	GAAGACATTGCCTGTGCTCGTG
Reverse Sequence:	GGTTTGCTGGTCTATCCTGTCG
ACCN:	<u>NM_001085430</u> , <u>NM_001085430.1</u> , <u>NM_001085430.2</u> , <u>BC041467</u> , <u>BC093905</u> , <u>BC093907</u> , <u>BU075516</u>
UniProt ID:	<u>Q0P5P2</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)