

## Product datasheet for **HP204296**

### **FAM162B Human qPCR Primer Pair (NM\_001085480)**

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221303
Forward Sequence:	CCGCAAGAAACAAAGCTCGAGTG
Reverse Sequence:	GCAGCTTCTTCACGCCACTTAG
ACCN:	<u>NM_001085480</u> , <u>NM_001085480.1</u> , <u>NM_001085480.2</u> , <u>BC038997</u> , <u>NM_001085480.3</u>
UniProt ID:	<u>Q5T6X4</u>
Synonyms:	bA86F4.2; C6orf189
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)