

## Product datasheet for **HP204300**

### TMEM150B Human qPCR Primer Pair (NM\_001085488)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284417
Forward Sequence:	GTCTGCACCATCCTCATTGTGG
Reverse Sequence:	CGGAGAAGTCAACGGCTAAGAG
ACCN:	<u><a href="#">NM_001085488</a></u> , <u><a href="#">NM_001085488.1</a></u> , <u><a href="#">NM_001085488.2</a></u> , <u><a href="#">BC020862</a></u> , <u><a href="#">BC171893</a></u> , <u><a href="#">BC171901</a></u> , <u><a href="#">BC171904</a></u> , <u><a href="#">NM_001085488.3</a></u>
UniProt ID:	<u><a href="#">A6NC51</a></u>
Synonyms:	DRAM3; TMEM224; TTN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)