

Product datasheet for **HP204332**

C5ORF44 (TRAPPC13) Human qPCR Primer Pair (NM_001093755)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80006
Forward Sequence:	GCATCCAAGTTCTTCGCTCTGTC
Reverse Sequence:	CCACACAGACTTGTGCGATGTC
ACCN:	<u>NM_001093755</u> , <u>NM_001093755.1</u> , <u>BC012006</u> , <u>BC047475</u> , <u>BC117128</u> , <u>BC142994</u> , <u>BC143691</u> , <u>BC143693</u> , <u>BE222610</u> , <u>BM664695</u> , <u>BX641157</u>
UniProt ID:	<u>A5PLN9</u>
Synonyms:	C5orf44; SHLD3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)