

Product datasheet for **HP204337**

FAM195B (MCRIP1) Human qPCR Primer Pair (NM_001093767)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	348262
Forward Sequence:	ACGAGGAGAACGTCCGCTTCAT
Reverse Sequence:	AGGGACCTTCTCCACATACTCC
ACCN:	<u>NM_001093767</u> , <u>NM_001093767.1</u> , <u>NM_001093767.2</u> , <u>BC111795</u> , <u>BC017108</u> , <u>BC019318</u> , <u>BC041068</u> , <u>BC063557</u> , <u>BC094876</u> , <u>BE831904</u> , <u>BG830426</u> , <u>BU188846</u> , <u>BX649156</u> , <u>NM_001093767.3</u>
UniProt ID:	<u>C9JLW8</u>
Synonyms:	FAM195B; GRAN2; MCRIP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)