

Product datasheet for **HP204338**

SFTPA1 Human qPCR Primer Pair (NM_001093770)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	653509
Forward Sequence:	CACCTGGAGAAATGCCATGTCC
Reverse Sequence:	AAGTCGTGGAGTGTGGCTTGA
ACCN:	<u>NM_001093770</u> , <u>NM_001093770.1</u> , <u>NM_001093770.2</u> , <u>BC111570</u> , <u>BC029913</u> , <u>BC171875</u> , <u>BI820937</u> , <u>BM996243</u> , <u>NM_001093770.3</u>
UniProt ID:	<u>Q8IWL2</u>
Synonyms:	COLEC4; PSAP; PSP-A; PSPA; SFTP1; SFTPA1B; SP-A; SP-A1; SP-A1 beta; SP-A1 delta; SP-A1 epsilon; SP-A1 gamma; SPA; SPA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)