

## Product datasheet for **HP204347**

### RPP38 Human qPCR Primer Pair (NM\_001097590)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10557
Forward Sequence:	GGACGAAGTAAGAGCCATCATCC
Reverse Sequence:	GGCTTTCCAGAGGTTCACTCTC
ACCN:	<u><a href="#">NM_001097590</a></u> , <u><a href="#">NM_001097590.1</a></u> , <u><a href="#">NM_001097590.2</a></u> , <u><a href="#">BC029494</a></u> , <u><a href="#">BC029494.1</a></u> , <u><a href="#">BP247818</a></u> , <u><a href="#">BU180901</a></u> , <u><a href="#">BU197354</a></u> , <u><a href="#">NM_001097590.3</a></u>
UniProt ID:	<u><a href="#">P78345</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)