

Product datasheet for **HP204358**

oncomodulin (OCM) Human qPCR Primer Pair (NM_001097622)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	654231
Forward Sequence:	CGGTTCATAGACAACGACCAGAG
Reverse Sequence:	CATTATCCGCCGAGCCATCAA
ACCN:	<u>NM_001097622</u> , <u>NM_001097622.1</u> , <u>BC113706</u> , <u>BC069468</u> , <u>BC069512</u> , <u>BC069550</u> , <u>BC113704</u> , <u>BC147028</u> , <u>BC147029</u> , <u>NM_001097622.2</u>
UniProt ID:	<u>P0CE72</u>
Synonyms:	OCM1; OM; ONCM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)