

Product datasheet for **HP204377**

GAGE12F Human qPCR Primer Pair (NM_001098405)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100008586
Forward Sequence:	GGAACCAGCAACTCAACGTCAG
Reverse Sequence:	GGACCATCTTCACACTCACACC
ACCN:	<u>NM_001098405</u> , <u>NM_001098405.1</u> , <u>NM_001098405.2</u>
UniProt ID:	<u>O76087</u>
Synonyms:	AL4; CT4.7; GAGE-7; GAGE-7B; GAGE-8; GAGE7; GAGE7B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)