

## Product datasheet for **HP204403**

### CEACAM21 Human qPCR Primer Pair (NM\_001098506)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	90273
Forward Sequence:	GGGAGTATCAGTGTGAGGTCTC
Reverse Sequence:	CCGATCAGGATGCCTAGAGTGT
ACCN:	<u><a href="#">BC012001</a></u> , <u><a href="#">NM_001098506</a></u> , <u><a href="#">NM_001098506.1</a></u> , <u><a href="#">NM_001098506.2</a></u> , <u><a href="#">NM_001098506.3</a></u> , <u><a href="#">BC106727</a></u> , <u><a href="#">BI833768</a></u> , <u><a href="#">BI913326</a></u> , <u><a href="#">NM_001098506.4</a></u>
UniProt ID:	<u><a href="#">Q3KPI0</a></u>
Synonyms:	CEACAM3; R29124_1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)