

Product datasheet for **HP204416**

CKAP2 Human qPCR Primer Pair (NM_001098525)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26586
Forward Sequence:	CCGTTGACCAGCGAAGACATAC
Reverse Sequence:	TTAGCACTCTTCCTTGCCAGCT
ACCN:	<u>NM_001098525</u> , <u>NM_001098525.1</u> , <u>NM_001098525.2</u> , <u>BC130296</u> , <u>BC010901</u> , <u>BC066920</u> , <u>BC105806</u> , <u>BC136332</u> , <u>BM541375</u>
UniProt ID:	<u>Q8WWK9</u>
Synonyms:	LB1; se20-10; TMAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)