

Product datasheet for **HP204441**

RBM47 Human qPCR Primer Pair (NM_001098634)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54502
Forward Sequence:	ATCCCTACCGTCAACCCAGTTG
Reverse Sequence:	GCTCCACTGTGTGGATTTTGCC
ACCN:	<u>NM_001098634</u> , <u>NM_001098634.1</u> , <u>BC126261</u> , <u>BC034402</u> , <u>BC071585</u> , <u>BC143942</u> , <u>BM837005</u> , <u>BP332501</u> , <u>BX648404</u>
UniProt ID:	<u>A0AV96</u>
Synonyms:	NET18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)