

Product datasheet for **HP204452**

GTF2H2C Human qPCR Primer Pair (NM_001098728)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	728340
Forward Sequence:	TCTTTAGCAGCCTTACAACCTTGCG
Reverse Sequence:	CAGTGCAAACGCGAACTTCTGC
ACCN:	NM_001098728 , NM_001098728.1 , BC105980 , BC157894 , BC171860 , NM_001098728.3
UniProt ID:	Q6P1K8
Synonyms:	GTF2H2C_2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)