

Product datasheet for **HP204490**

ZNF862 Human qPCR Primer Pair (NM_001099220)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 643641 |
| Forward Sequence: | TCACTGATGGCATCCACAGCTC |
| Reverse Sequence: | CCTCAAACACCACAGGAACCTC |
| ACCN: | <u>NM_001099220</u> , <u>NM_001099220.1</u> , <u>NM_001099220.2</u> , <u>BC041427</u> , <u>BC045757</u> , <u>BC092506</u> , <u>BC113087</u> , <u>BQ270712</u> , <u>BU164247</u> , <u>NM_001099220.3</u> |
| UniProt ID: | <u>O60290</u> |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)