

## Product datasheet for **HP204545**

### RPL22L1 Human qPCR Primer Pair (NM\_001099645)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	200916
Forward Sequence:	CAAGTGGATGGACCCCTTCTCC
Reverse Sequence:	GGCACATCATTCTTGCTCTCC
ACCN:	<a href="#">NM_001099645</a> , <a href="#">NM_001099645.1</a> , <a href="#">BC107708</a> , <a href="#">BC062731</a> , <a href="#">BM451642</a> , <a href="#">BU540139</a> , <a href="#">NM_001099645.2</a>
UniProt ID:	<a href="#">Q6P5R6</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)