

## Product datasheet for **HP204755**

### FMN1 Human qPCR Primer Pair (NM\_001103184)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 342184  |
| Forward Sequence: | CATGAGAAGGAGTCGCTAAGAGC   |
| Reverse Sequence: | TGAAGTCTGCCAGGAGTCCTGT  |
| ACCN:             | <u><a href="#">NM_001103184</a></u> , <u><a href="#">NM_001103184.1</a></u> , <u><a href="#">NM_001103184.2</a></u> , <u><a href="#">NM_001103184.3</a></u> , <u><a href="#">BC029107</a></u> , <u><a href="#">BC101971</a></u> , <u><a href="#">BC103692</a></u> , <u><a href="#">BC107593</a></u> , <u><a href="#">BC172380</a></u> , <u><a href="#">BM663620</a></u> , <u><a href="#">NM_001103184.4</a></u> |
| UniProt ID:       | <u><a href="#">Q68DA7</a></u>   |
| Synonyms:         | FMN; LD   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)