

Product datasheet for **HP204839**

ADAM8 Human qPCR Primer Pair (NM_001109)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	101
Forward Sequence:	TGCTGGAGGTGGTGAATCACGT
Reverse Sequence:	TCAGGAGTTCTCCAGTGTGAC
ACCN:	<u>NM_001109</u> , <u>NM_001109.1</u> , <u>NM_001109.2</u> , <u>NM_001109.3</u> , <u>NM_001109.4</u> , <u>BC064500</u> , <u>BC115404</u> , <u>BI908229</u> , <u>BM144968</u> , <u>BM146234</u> , <u>BQ007791</u> , <u>NM_001109.5</u>
UniProt ID:	<u>P78325</u>
Synonyms:	CD156; CD156a; MS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)