

## Product datasheet for **HP204844**

### PTPRB Human qPCR Primer Pair (NM\_001109754)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5787
Forward Sequence:	TCTTCCCGACAAGTGGTTGTGG
Reverse Sequence:	AGCCAGGAAACGCTGAGGTAGT
ACCN:	<u><a href="#">NM_001109754</a></u> , <u><a href="#">NM_001109754.1</a></u> , <u><a href="#">NM_001109754.2</a></u> , <u><a href="#">NM_001109754.3</a></u> , <u><a href="#">BC043182</a></u> , <u><a href="#">BC051329</a></u> , <u><a href="#">BC101679</a></u> , <u><a href="#">BC113463</a></u> , <u><a href="#">BC143356</a></u> , <u><a href="#">BC143360</a></u> , <u><a href="#">BE042873</a></u> , <u><a href="#">BQ710210</a></u> , <u><a href="#">BX647238</a></u> , <u><a href="#">BX648245</a></u> , <u><a href="#">BX648771</a></u> , <u><a href="#">NM_001109754.4</a></u>
UniProt ID:	<u><a href="#">P23467</a></u>
Synonyms:	HPTP-BETA; HPTPB; PTPB; R-PTP-BETA; VEPTP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)