

Product datasheet for **HP204845**

FP100 (FAAP100) Human qPCR Primer Pair (NM_001109760)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80233
Forward Sequence:	AGCCAGGATGACCACAGAGAGT
Reverse Sequence:	TTGTCAGTGCCTTGTTCCGCTG
ACCN:	<u>NM_001109760</u> , <u>NM_001109760.1</u> , <u>BC117139</u> , <u>BC001724</u> , <u>BC008883</u> , <u>BC021968</u> , <u>BC117141</u> , <u>BC143861</u>
UniProt ID:	<u>Q0VG06</u>
Synonyms:	chromosome 17 open reading frame 70; FAAP100; Fanconi anemia associated protein 100 kDa subunit; Fanconi anemia core complex 100 kDa subunit; FLJ22175; FLJ30151
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)