

Product datasheet for **HP204860**

EBF4 Human qPCR Primer Pair (NM_001110514)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57593
Forward Sequence:	TTGCAGGTCGTGTTCCGAAACG
Reverse Sequence:	CTGCTTGGAAGTTAGGAGAGG
ACCN:	<u>BC054347</u> , <u>NM_001110514</u> , <u>NM_001110514.1</u> , <u>BC054347.1</u> , <u>BC010557</u> , <u>BC019106</u> , <u>BC037555</u> , <u>BG750986</u>
UniProt ID:	<u>Q9BQW3</u>
Synonyms:	COE4; O/E-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)