

Product datasheet for **HP204945**

C17orf85 (NCBP3) Human qPCR Primer Pair (NM_001141118)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55421
Forward Sequence:	GGAACCTGGAGCAGCCAGAAGA
Reverse Sequence:	CCGCTGAATACGACGGGAATGA
ACCN:	<u>NM_001141118</u> , <u>NM_001141118.1</u> , <u>NM_001141118.2</u> , <u>BC010707</u> , <u>BC011733</u> , <u>BC146419</u> , <u>BC148856</u> , <u>BM974475</u> , <u>BP344033</u> , <u>NM_001141118.3</u>
UniProt ID:	<u>Q53F19</u>
Synonyms:	C17orf85; ELG; HSA277841
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)