

Product datasheet for **HP204971**

JMJD7 Human qPCR Primer Pair (NM_001114632)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100137047
Forward Sequence:	GGAGTCCTCTATGTGCAGAAGC
Reverse Sequence:	CAGCCAGAAGTTCACAGCATCG
ACCN:	BC025290 , NM_001114632 , NM_001114632.1 , BC025290.1 , BC067834 , BC110415 , BC139899 , BU739980 , NM_001114632.2
UniProt ID:	P0C870
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)