

Product datasheet for **HP204975**

PLCL1 Human qPCR Primer Pair (NM_001114661)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5334
Forward Sequence:	CGGTTCATCCATTACGAGAAGCC
Reverse Sequence:	AACAGGCACTGCTTCAGACTGG
ACCN:	<u>NM_001114661</u> , <u>NM_001114661.1</u> , <u>BC111985</u> , <u>BC101531</u> , <u>BP340710</u> , <u>BX537442</u>
UniProt ID:	<u>Q15111</u>
Synonyms:	MGC126580; MGC138190; OTTHUMP00000163640; phospholipase C, epsilon; phospholipase C-like 1; PLC-L; PLCE; PLCL; PLDL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)